

Laser Tissue Blood Flowmeter

OMEGAFLO Serise




Non-contact FLO-N1
 Blood flow measurement without contacting the Probe
 Blood flow value is averaged over wider area. Good reproducibility.

Contact FLO-C1
 Conventional laser tissue blood flowmeter
 Many kinds of probes are available.
 Specialized probes can be ordered.

Continuously and non-invasively

Measurement without contacting the probe



FLO-N1

[Outline]

OMEGAFLO measures tissue blood flow, blood volume and flow velocity continuously and non-invasively. It is used for evaluation of the effect of drugs, transplantations and physiological research.



FLO-N1, the non-contact type, can measure- tissue blood flow without contacting the probe on tissue, and it allows measurement and operations on the tissue at the same time. The measurement depth- is about 1 mm from the tissue surface.

[Theory]


The tissue under study is irradiated by the laser light through the optical fiber probe, and the light is scattered many times in the tissue. Some of the photons are scattered by erythrocytes in the tissue, and light is shifted in frequency related to the moving velocity of the erythrocytes.

Also , the amplitude of the scattered light is proportional to the volume of the erythrocytes in the tissue. The blood flow signals are obtained from the frequency shift and its amplitude.

[Type ST-N probe]



[Type DS probe]



Specification			
Model		FLO-N1	FLO-C1
Method		Non-contact	Contact
Laser Power at Probe		< 3 mW	< 2 mW
Measurement Area and Depth		Inside of the 15mm x 1mm	about 1mm x 1mm
laser		780nm	
Parameters		Tissue Blood Flow(FLOW) Tissue Blood Volume(MASS) Tissue Blood Velocity(VEL.)	
Indicator		FLOW, MASS, VELOCITY (Received Power)	
Output		FLOW, MASS, VELOCITY 0-10V Analog	
Time Constant		0.1/3.0/5.0 Sec	
Input Voltage		A C 100-220V	
Size		256(W)_324(D)_69mm	
Weight		4 k g	



<http://www.omegawave.co.jp>

Head Office : Sun noble 2-20-3 katamachi, Fuchu, Tokyo, Japan